

**MBSIREPI\_ref**

TA: 5.0 s PM: REF Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	20
Dist. factor	25 %
Position	L0.0 A3.6 H30.3 mm
Orientation	T > C-5.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-7

**Contrast - Common**

TR	50.00 ms
TE	27.60 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

**Resolution - Common**

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm

**Resolution - Common**

Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	20
Dist. factor	25 %
Position	L0.0 A3.6 H30.3 mm
Orientation	T > C-5.5
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 A3.6 H30.3 mm
Orientation	T > C-5.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A3.6 H30.3
L	0.0 mm
A	3.6 mm
H	30.3 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-5.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H

## System - Miscellaneous

Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

## System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

## System - Adjust Volume

Position	L0.0 A3.6 H30.3 mm
Orientation	T > C-5.5
Rotation	0.00 deg
A >> P	210 mm
R >> L	210 mm
F >> H	99 mm
Reset	Off

## System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

## System - Tx/Rx

Frequency 1H	123.195810 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

## Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	1

## Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

## Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

## Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

## Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

## Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

## Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

## Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher time [1]	0.0 ms,ms
crusher time [2]	0.4 ms,ms
crusher time [3]	22.00 ms,ms

## Sequence - Assistant

Mode	Off
------	-----